



*Invacare*<sup>®</sup>

# ***Essential Basic***


Cost effective comfort mattress



The **Essential Basic** mattress offers the client a comfortable and supportive surface. This is an ideal entry level mattress for individuals who do not require therapy for pressure ulcer prevention.

- **Multi-stretch vapor-permeable, water resistant polyurethane fabric cover.**
- **Zipped around 3 sides for ease of inspection and cleaning.**
- **Meets fire safety standards.**

## Features and Options



### Entry Level



Ideal for clients who do not require therapy for pressure ulcer prevention.

### Cover

Multi-stretch polyurethane fabric cover.


## Technical data

For more information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Size range	Max. user weight	Warranty
			1 year
<b>Essential Basic</b>	2000 x 880 x 150 mm 2050 x 880 x 150 mm	140 kg (22 stone)	

## Colour

Blue



Mattresses is manufactured to fire safety standard:  
BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

